



# Viral Transport Media (VTM) with Swab

**Catalog # 52018 (Includes 52010 and 52014)**

- Ready-to-use, compatible transport media to safely and conveniently meet all of your viral testing needs — including SARS-CoV-2
- Manufactured according to FDA requirements under KSL FDA establishment license and cGMP compliant quality management system
- Made in the USA — transport medium tubes are manufactured and packaged by KSL in New York State and NYS-manufactured swabs are available
- Ready to ship today

**Detailed specifications on back.**



## Specifications

<b>Catalog Number</b>	<ul style="list-style-type: none"><li>• 52018 (Includes 52010 and 52014)</li></ul>
<b>Components</b>	<ul style="list-style-type: none"><li>• Sterile VTM Filled Tube</li><li>• Nylon Flocked Swab</li><li>• Packaged inside a sealed pouch</li></ul>
<b>Tube</b>	<ul style="list-style-type: none"><li>• Volume: 10mL PP tube with an HDPE cap (the volume of the media is 3mL)</li><li>• Length – Body: 101mm</li><li>• Length – Body with Cap: 105mm</li><li>• Diameter: 15mm</li><li>• Sterile with a leak-proof cap</li></ul>
<b>Media</b>	<ul style="list-style-type: none"><li>• CDC Formulated Viral Transport Media</li><li>• Sterile</li></ul>
<b>Swab</b>	<ul style="list-style-type: none"><li>• Sterile flocked nylon nasopharyngeal/oropharyngeal swab</li><li>• Packaged inside a sealed pouch</li></ul>
<b>Packaging</b>	<ul style="list-style-type: none"><li>• Each unit consists of a media-filled tube with a swab packaged inside a sealed pouch</li><li>• Pouches are sold in bags of 50 units</li></ul>
<b>Labeling</b>	<ul style="list-style-type: none"><li>• Vinyl-impregnated water/temperature resistant 2"x1" labels</li><li>• FDA-compliant UDI labels</li></ul>
<b>Stability</b>	<ul style="list-style-type: none"><li>• One year from date of manufacture when stored between 2° and 8°C</li><li>• Six months when stored between 2° and 30°C</li></ul>
<b>Storage</b>	<ul style="list-style-type: none"><li>• Between 2° and 8°C or between 2° and 30°C (see "Stability" above)</li><li>• Do not freeze</li></ul>



## To Order, Contact KSL Biomedical:

1000 Youngs Road  
Suite 210  
Buffalo NY 14221 USA

1.800.960.1080  
orders@kslbiomedical.com  
www.kslbiomedical.com

