

# MONO Rapid Test Card

## (Whole Blood/Serum/Plasma)



**Infectious Mononucleosis (IM)** is caused by the Epstein-Barr virus, which is a member of the herpesvirus family. Symptoms of IM are fever, sore throat and swollen lymph glands. In very rare cases, heart or central nervous system problems may occur. Diagnosis of IM is made based on the presence of heterophile antibodies. Infectious Mononucleosis heterophile antibodies belong to the IgM class. They are present in 80-90% of acute IM cases and can be detected in 60-70% of patients during the first week of clinical illness.



The **MONO Rapid Test Card** is a rapid chromatographic immunoassay for the qualitative detection of Infectious Mononucleosis heterophile antibodies in whole blood, serum or plasma as an aid in the diagnosis of Infectious Mononucleosis.

**Reading Time:** 5 Minutes

**Sensitivity:** 97.6%

**Specificity:** 97.8%

**Accuracy:** 97.7%

**Product Advantages:**

- Fast result
- Easy visual interpretation
- Simple operation
- High accuracy

Catalog #	Description
60102-50	MONO Rapid Test Card (whole blood/serum/plasma) – 50 Tests
60102-100	MONO Rapid Test Card (whole blood/serum/plasma) – 100 Tests

5100 South Service Road, Unit 18  
Burlington, Ontario, Canada L7L 6A5  
Tel (905) 333-8188 | Fax (905) 333-0500  
[order@pulsescientific.com](mailto:order@pulsescientific.com)